Revision Date: May 9, 2025

Cerebral Corridor Creator Instructions for Use

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Content

1. Device Description	3
2. Device Specifications and Size Parameters	4
3. Selection of Cerebral Corridor Creator	4
4. Intended Use	5
5. Indications	5
6. Contraindications	5
7. Precautions and Warning	5
8. Directions for Use	6
9. Expected clinical benefit	9
10. Shelf Life	10
11. Storage, Transportation and Disposal	10
12. Package and Labeling	10
13. Links to the Summary of Safety and Clinical Performance	10
14. Disclaimer of Warranty and Limitation of Remedy	10

1. Device Description

The Cerebral Corridor Creator consists of a Cerebral Corridor Creator tube and a Cerebral Corridor Creator balloon.

- **Description of Cerebral Corridor Creator**


The Cerebral Corridor Creator (CCC) is a specially designed medical device for neurosurgeons. It provides a way to access surgical targets in brain without unduly damaging surrounding tissues. The CCC consists of a balloon and a tube. A tubular corridor is created and the collateral trauma can be minimized along with the balloon inflated slowly and repetitively by a syringe. And the tube can make the corridor stability and reliability, which can also protect surrounding tissues during the operation. With the goal of establishing a working corridor in minor brain intervention, the CCC can be used in neurosurgery such as hematoma, hydrocephalus, tumor and other lesions etc.

- **Product Structure**

a) Balloon: is intended for brain tissue retraction and corridor creating.

The balloon consists of an inner tube, an outer tube, a latex balloon, a three-way stopcock and a connector. (see Fig. 1)

b) Tube: is intended for maintaining the corridor. (see Fig. 2)

1. Latex balloon 2.Outter tube 3. Three-way stopcock 4.Inner tube 5.Connector

Fig.1 Structure of Balloon


Fig. 2 Structure of Tube

2. Device Specifications and Size Parameters

2.1. The specifications and size parameters of Cerebral Corridor Creator are shown in Table 1.

Table 1 Specifications of Cerebral Corridor Creator

Part Name	Model	O.D.	Working Diameter (mm ± 2 mm)	Overall Length	Working Length (mm ± 5 mm)	Recomm. Air Volume (mL)
Cerebral Corridor Creator Balloon	QK-18L22	8F (≤2.6mm)	18	22±1cm	50	10
	QK-17L22	8F (≤2.6mm)	17		50	9
	QK-16L22	8F (≤2.6mm)	16		50	8
	QK-15L22	8F (≤2.6mm)	15		50	7
	QK-14L22	8F (≤2.6mm)	14		50	6
	QK-13L22	8F (≤2.6mm)	13		50	5
	QK-12L22	8F (≤2.6mm)	12		50	4
	QK-10L24	8F (≤2.6mm)	10	24±1cm	50	4
	QK-08L24	8F (≤2.6mm)	8		50	3
	QK-06L24	8F (≤2.6mm)	6		50	2
Cerebral Corridor Creator Tube	DG-5015B	33±2mm	/	50±5mm	/	
	DG-7013B			70±5mm		

2.2. Product performance

- Visual Inspection: The Cerebral Corridor Creator Balloon and Cerebral Corridor Creator Tube should be clear as no particular and foreign matter. The Cerebral Corridor Creator Balloon and Cerebral Corridor Creator Tube should be no obvious defects as burr, deformation, break, etc.
- Freedom From leakage: The catheter shall not leak.
- EO residuals: EO Sterilization is adopted. The EO residuals should be less than 10ug/g.

3. Selection of Cerebral Corridor Creator

Firstly, the length of the tube can be determined on the basis of depth of lesion which can be estimated by preoperative imaging examination, while the inner diameter of the tube can be determined on the basis of tubular corridor's size which can be estimated on the lesions size combination with surgeon's operating skills and proficiency, after that, the corresponding balloon specification can be determined according to Table 1 the Specifications of Cerebral Corridor Creator in section 2 in IFU.

Considering of safety, the tubular corridor's size should not be more than 22 mm. If it exceeds 22mm, brain tissue may be damaged.

4. Intended Use

Cerebral Corridor Creator is intended for use in the neurosurgery, as a brain tissue retractor and holder device, provides a working corridor to neurosurgery operation.

5. Indications

The Cerebral Corridor Creator is intended to be used for operating on patients for the following situations:

- Hypertensive cerebral hemorrhage (basal ganglia, lobes, thalamus, ventricles).
- Biopsy of patients with obstructive hydrocephalus.
- Retractor with various auxiliary operation methods..

6. Contraindications

- Cerebral hemorrhage caused by intracranial aneurysm deformity, intracranial aneurysm, hemopathy, etc
- Vascular malformations and aneurysm of uncertain size, shape and location;
- Abundant vascular tumor.
- Cerebral hernia.
- Coagulation disorders.
- Cerebral traumatic hemorrhage
- Other cases not suitable for the device considered by surgeon.

7. Precautions and Warning

- It will cause balloon leakage or rupture if there is sharp object inside the working channel. Immediately withdraw and replace the balloon if the abnormal situations happen to the balloon when inflation.
- This product contains natural rubber latex, which may cause allergic reactions such as: itching, redness, urticaria, swelling, fever, dyspnea, asthma-like symptoms, hypotension, shock, etc. If one or more of these symptoms are observed in the patient, stop procedure immediately and take appropriate treatment.

- You should read the instructions before use and make sure all parts are not damaged. Don't use the product if its package is damaged. Use within its valid period.
- DO NOT RESTERILIZE AND/OR REUSE. After being used, dispose of the devices in accordance with hospital, administrative, and/or government policy. Reuse or resterilization may create a risk of contamination of the device and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient. Reuse or resterilization may compromise the structural integrity of the device and/or lead to device failure, which, in turn, may result in patient injury, illness, or death. SYM company will not be responsible for any direct, incidental, or consequential damages resulting from resterilization or reuse.
- The product should be used by the trained and qualified doctors; Doctors should be strictly screened for indications.
- Suggestion that surgery can be only performed in hospitals which is qualified for surgery at any time. The doctor must be trained qualified in professional microscope or endoscopic surgery, and master the technology of the Department of Neurosurgery, and complete 5 cases of surgery under guidance of the technical maturity doctors.
- The packaged product should be stored in the clean room which is under no more than 80% RH, room temperature, cool, dry, dark, no corrosive gas and ventilated.
- Carefully check the product before use, prepare backup spare parts for abnormal situations.
- The user and/or patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State.

8. Directions for Use

The recommended procedures are shown as following:

8.1. Preoperative Preparation

8.1.1. Perform intracranial vascular examination using techniques such as MRA, DSA, CTA, and etc to eliminate contraindications and determine lesion location and also the specification of cerebral corridor creator.

8.1.2. Perform a general anesthetic to patients.

8.1.3. Perform brain needle puncture at point of determined lesion location after craniotomy.

8.1.4. Unpack the package of cerebral corridor creator after check its integrity.

Warning: DO NOT use the device if its sterile package is open.

8.1.5. Flushing

Fully flushing the Cerebral Corridor Creator including balloon, inner tube and outer tube completely with normal saline (10-100u/ml) before using.

8.1.6. Exhaust air from the balloon with normal saline

Connect the syringe (we suggest rotary injection syringe) to the connector. Turn the valve of connector on. Inject normal saline thru the connector into the balloon. Make sure the balloon is placed at a low level while injecting 0.9% saline slowly. Deflate the balloon when it has been dilated. Repeat this procedure until there is no air remaining in the balloon lumen.

Caution: If there is any remaining air in the balloon lumen, the balloon may fail to expand to its required diameter. If balloon ruptures unexpectedly, it might cause some brain tissue damage. To avoid this, make sure to exhaust air from the balloon lumen completely.

8.1.7. Extending balloon

Push the inner tube straight forward until the other end of the inner tube is locked in its position. In this way the balloon is fully extended (see Fig.3).


Fig.3 Extending balloon

Caution: Don't turn inner tube in any direction while you are pushing forward the inner tube. That may make balloon twist around itself.

8.2. Procedures in operation

8.2.1. Insert extended balloon into predetermined position

First Connect balloon to the syringe and make sure the volume of normal saline does not exceed the balloon's maximal capacity. The balloon's maximal capacity is shown in Table 2. Then insert the balloon into the brain tissue according to the preoperative guided route.

Table 2. The balloon's maximal capacity

Specification	Maximal capacity	The nominal capacity of the syringe is recommended	Specification	Maximal capacity	The nominal capacity of the syringe is recommended
QK-06L24	2mL	2mL	QK-14L22	6mL	10mL
QK-08L24	3mL	5mL	QK-15L22	7mL	10mL
QK-10L24	4mL	5mL	QK-16L22	8mL	10mL
QK-12L22	4mL	5mL	QK-17L22	9mL	10mL
QK-13L22	5mL	5mL	QK-18L22	10mL	10mL

Caution: Make sure the volume of normal saline in the syringe does not exceed the balloon's maximal capacity.

8.2.2. Expanding surgical route with balloon to form a tubular corridor

Inject and withdraw normal saline into balloon repeatedly to inflate and deflate the balloon cautiously and slowly, and insert the balloon step by step into the preoperative guided route which punctured by the brain needle, in this way a tubular corridor can be created at the predetermined position. If any unexpected incident occurs, stop immediately and take back up methods.

8.2.3. Inserting the tube together with the balloon

After you create a tubular corridor, place the tube outside of the distal part of balloon. Inflate the balloon and then insert both tube and balloon slowly creep into the tubular corridor (Fig.4). Inflating and deflating the balloon appropriately would be helpful.


Fig. 4 Inserting the tube together with the balloon

8.2.4. Withdrawing

Withdraw the normal saline and deflate the balloon after the tube reaches the predetermined position. Pull the tube out slowly when the surgery is finished.

8.3. Post-procedure

- Re-examination is needed within 24 hours post operation.
- Antibiotic therapy.

9. Expected clinical benefit

According to clinical trial results, the Cerebral Corridor Creator can reduce intraoperative blood loss, and no adverse events occurred within 24 hours after surgery.

10. Shelf Life

The Cerebral Corridor Creator is sterilized by ethylene oxide. Please use the device before the “USE by” date specified on the label.

11. Storage, Transportation and Disposal

- The product should be kept in protective, well ventilated and clean room with no aggressive gas. Stored at room temperature. The upper limit for relative humidity is 80%.
- The product should be placed upward, and could not mix with heavy, corrosive, sharp objects, and prevent from direct sunlight and moisture.
- While in transportation, the product should be carefully handled, and prevent from being squeezed and exposure under direct sunlight and moisture.
- The product is disposable. Please contact with SYM or distributor directly if you have questions to the product before use.
- The expired date is indicated on the product label. Do not use the product if it is expired.

12. Package and Labeling

The Cerebral Corridor Creator is supplied sterile. The product package consists of a sterilized pack, an outer paper box and a transportation package. All the components are sealed in sterilized pack on which a label and a sterilization indicator are stuck. The product is sterilized and then put in the outer box made of art paper with IFU. During transportation, the product with outer box will be packaged in a transportation package made of corrugated paper to protect the product from damage during transportation.

13. Links to the Summary of Safety and Clinical Performance

The SSCP will be available via Eudamed and can be located by using the Basic UDI of the Cerebral Corridor Creator. The Basic UDI for Cerebral Corridor Creator is 69396986SYM0001EN2D.

14. Disclaimer of Warranty and Limitation of Remedy

- SYM is not responsible for any damage caused due to improper handling, operation and storage, including use after the expiry date printed on the label.
- SYM is not responsible for any damage to the disposable Cerebral Corridor Creator caused by transportation, handling and storage at the health care facility, whether physical damage or human injury.
- SYM is not responsible for any damage caused by reprocessing or reuse of the disposable Cerebral Corridor Creator, whether physical damage or human injury.
- Do not use the disposable Cerebral Corridor Creator with other manufacturer's products. SYM is not responsible for any impairment to structural integrity and/or function which may arise from the use of the disposable Cerebral Corridor Creator with other manufacturer's products.
- SYM is not responsible for any damage to patients and users caused by unintended use of the disposable Cerebral Corridor Creator.

Symbol	Symbol Definitions
	Manufacturer
	Authorized representative
	Country of manufacture
	Model number
	Use-by date
	Batch code
	Serial number
	Sterilized using ethylene oxide
	Do not resterilize
	Do not use if package is damaged and consult instructions for use

	Do not reuse
	Consult instructions for use or consult electronic instructions for use
	Caution
	Contains or presence of natural rubber latex
	Medical device
	Unique device identifier
	Single sterile barrier system
	Single sterile barrier system with protective packaging inside
	Fragile, handle with care
	Keep away from sunlight
	Keep dry
	Usable length of the device
	CE marking 2797: Identification number of notified body

Manufacturer:

Shenzhen Shineyard Medical Device Company Limited

3/F, Changfeng Industrial Block, No.3 Liuxian Road, Xin'an, Bao'an District,

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SUNGO Europe B.V.